

## **Assay method for Protease activity (Folin method pH 8.0)**

### **Solutions**

1) 0.2 mol/L Hydrochloric acid solution

Dissolve 18 mL of Hydrochloric acid in about 500 mL of water and dilute to 1000 mL with water.

2) 6 mol/L Acetic acid solution

Dilute 180.0 g of Acetic acid to 500 mL with water.

3) Sodium hydroxide solution

Dissolve 43.0 g of Sodium hydroxide in about 400 mL of water and cool. Dilute to 1000 mL with water.

4) 0.55 mol/L Sodium carbonate solution

Dissolve 116.58 g of Anhydrous sodium carbonate in about 1500 mL of water and dilute to 2000 mL with water.

5) 1 mg/mL Tyrosine stock solution

Dissolve 0.1000 g of L-Tyrosine in 0.2 mol/L Hydrochloric acid solution and dilute to 100 mL with 0.2 mol/L Hydrochloric acid solution.

6) Trichloroacetic acid solution

Dissolve 36.00 g of Trichloroacetic acid and Anhydrous sodium acetate in about 1600 mL of water. Add to 110 mL of 6 mol/L Acetic acid solution and dilute to 2000 mL with water.

7) 0.1 mol/L Potassium dihydrogen phosphate solution

Dissolve 27.21 g of Potassium dihydrogen phosphate in about 1000 mL of water and dilute to 2000 mL with water.

8) 0.1 mol/L Disodium hydrogen phosphate solution

Dissolve 28.39 g of Anhydrous disodium hydrogen phosphate about 1000 mL of water and dilute to 2000 mL with water.

9) 0.1 mol/L Phosphate buffer pH 8.0

Add 0.1 mol/L Potassium dihydrogen phosphate solution to 0.1 mol/L Disodium hydrogen phosphate solution until the pH stabilizes at 8.00

10) 0.02 mol/L Phosphate buffer pH 8.0

Mix 2 L of 0.1 mol/L Phosphate buffer pH 8.0 and 8 L of water.

11) 3 fold Folin & Ciocalteu's phenol reagent

Dilute 50 mL of Folin & Ciocalteu's phenol reagent to 150 mL with water.

12) 0.05 mol/L Disodium hydrogen phosphate solution

Dissolve 14.2 g of Anhydrous disodium hydrogen phosphate in about 1000 mL of water and dilute to 2000 mL with water.

13) Substrate solution

Dissolve 1.20 g of Hammarsten casein in 160 mL of 0.05 mol/L Disodium hydrogen phosphate solution while heating at 60 °C and stirring occasionally. After cooling, add Sodium hydroxide solution until the pH stabilizes at 8.00. Dilute to 200 mL with water.

### Preparation of the sample solution

Dissolve the sample in an appropriate amount of water or a buffer solution.

## PROCEDURE

### Tyrosine Standard Curve

Dilute 1 mL, 2 mL, 3 mL and 4 mL of 1 mg/mL Tyrosine stock solution to 100 mL with 0.2 mol/L Hydrochloric acid solution. Pipet 2 mL of each these solutions in a test tube. Add 5 mL of 0.55 mol/L Sodium carbonate solution and 1 mL of 3 fold Folin&Ciocalteu's phenol reagent, and mix. Place the tubes in water bath maintained at  $37 \pm 0.5$  °C for 30 minutes. Measure the absorbance at 660 nm.

As the blank, use 0.2 mol/L Hydrochloric acid solution instead of diluted Tyrosine stock solution in the same manner.

Plot the absorbance against the quantity of Tyrosine ( $\mu$ g) to determine the slop of the standard curve.

$$F = \{20/(A_{s20}-A_{sBL})+40/(A_{s40}-A_{sBL})+60/(A_{s60}-A_{sBL})+80/(A_{s80}-A_{sBL})\}/4$$

### Sample Assay

Pipet 5 mL of Substrate solution into test tubes and place in a water bath maintained at  $37 \pm 0.5$  °C for 10 minutes. Add 1mL of sample solution and mix. Place the tubes in the water bath at  $37 \pm 0.5$  °C for 10 minutes. Add 5mL of Trichloroacetic acid solution and mix. Place the tubes in the water bath at  $37 \pm 0.5$  °C for 30 minutes and filtrate with filter paper (Toyo: No.131).

Pipet 5 mL of 0.55mol/L Sodium carbonate solution in test tubes, add 2 mL of above filtrate and 1 mL of 3 fold Folin&Ciocalteu's phenol reagent. Mix and place in the water bath at  $37 \pm 0.5$  °C for 30 minutes. Measure the absorbance at 660nm.

As the blank, pipet 1 mL of Sample solution in test tubes, add 5 mL of Trichloroacetic acid solution and mix. Add 5 mL of Substrate solution and mix. Place the tubes in the water bath at  $37 \pm 0.5$  °C for 30 minutes and filtrate with filter paper (Toyo: No.131).

Pipet 5m L of 0.55 mol/L Sodium carbonate solution in test tubes, add 2 mL of above filtrate and 1 mL of 3 fold Folin&Ciocalteu's phenol reagent. Mix and place in the water bath at  $37 \pm 0.5$  °C for 30 minutes. Measure the absorbance at 660 nm.

### Definition of Activity Unit

One unit is defined as the quantity of enzyme that liberates the coloring substance equivalent of 1 µg of Tyrosine for 1 minute under the conditions of the assay.

### CALCULATION

$$\text{Protease activity unit/g,mL} = (A_{10} - A_{0\text{BL}}) \times F \times (11/2) \times (1/10) \times n$$

$A_{10}$	: Absorbance of the Sample.
$A_{0\text{BL}}$	: Absorbance of the Blank.
F	: Quantity of Tyrosine from the calibration curve (µg)
11/2	: Conversion factor for the total volume
1/10	: Conversion factor for 1 minute
n	: Dilution factor of the enzyme